optilube.™

REF 1125

REF 1126

optilube.™

Instructions for use

optimummedical®

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UK Responsible Person and Distributor:

Optimum Medical Solutions Ltd
Tennant Hall, Blenheim Grove,
Leeds, LS2 9ET, United Kingdom
T: +44(0) 845 643 5479
enquiries@optimummedical.co.uk
www.optimummedical.co.uk

EU Importer:

Optimum Medical Solutions GmbH,
c/o Warth & Klein Grant Thornton AG,
Ulmenstraße 37-39,
60325 Frankfurt am Main, Germany
Email: enquiries@optimummedical.de

1125-1126-A0139-V1.0

STERILE R

MD

CE 2292

GBR INSTRUCTIONS FOR USE**PRODUCT REF/CONTENTS**1125 **OptiLube** 6ml Prefilled Syringe1126 **OptiLube** 11ml Prefilled SyringeThis is an **OptiLube** product information leaflet for the user.**OptiLube** has been prescribed for you by a Healthcare Professional. Do not share it with anyone else.

Read all of the product information leaflet prior to use and retain for your information.

If you are unsure about how to use this product, or whether it is suitable for you, please seek medical advice.

PRODUCT DESCRIPTION**OptiLube** is a sterile, single-use, pre-filled syringe with water-soluble jelly available in 6ml and 11ml syringes.**INGREDIENTS**100g of **OptiLube** contains:

- Purified Water;
- Propylene Glycol;
- Hydroxyethylcellulose.

INTENDED PURPOSE**OptiLube** is a medical device that has a lubricating effect that helps prevent trauma being caused to the patient during catheterisation procedures or other urethral procedures by effective lubrication.**OptiLube** is a clear, water-soluble, petroleum free lubricating jelly.**INDICATION FOR USE****OptiLube** is designed for lubrication of the urethra prior to insertion of urethral catheters and other urological medical devices including cystoscopes. It is also suitable for supra pubic catheterisation.**CONTRAINDICATIONS****OptiLube** must not be used in patients:

- with known hypersensitivity to any of the ingredients.

Do not use for intravenous (IV) or intramuscular (IM) injections.

Do not use orally. If swallowed, please seek medical advice.

Do not use in eyes. If it comes into contact with the eyes, please seek medical advice.

WARNINGS AND PRECAUTIONS**OptiLube** must only be used under the supervision of Healthcare Professionals in accordance with local guidelines, policies and procedures.

Keep out of reach of children.

Ensure **OptiLube** is suitable for the intended use and compatible with other medical devices to be used in conjunction with it.

This is a single-use device. Re-use of this device may result in patient infection/ cross-contamination.

Re/sterilisation, reprocessing, cleaning and disinfection may also compromise the product characteristics, resulting in trauma or infection to the patient.

OptiLube is latex free.**UNDESIRABLE SIDE EFFECTS**In rare cases local hypersensitivity reactions may occur, such as redness, stinging, blistering or itching. If this occurs, stop using **OptiLube** and consult your Healthcare Professional.Any serious incident or malfunction which may affect the safety of **OptiLube** should be reported to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.**HOW TO USE**The decision as to which size (6ml or 11ml) of **OptiLube** sterile lubricating jelly pre-filled syringe to be used is taken by the Healthcare Professional. Consult your Healthcare Professional for advice before using **OptiLube**. It is recommended that this product is used as part of an aseptic technique.

- Cleanse the opening to the urethra and surrounding area prior to use;
- Evenly peel back the paper backing and remove the sterile syringe;
- Remove the cap from the end of syringe;
- Apply a drop of jelly to urethral opening to make initial insertion easier;
- Insert the nozzle into the urethral opening and press the plunger of the syringe slowly to release the necessary amount of **OptiLube**. NOTE: **OptiLube** should be applied to the urethra, not directly onto the device;
- The lubrication characteristics of the jelly start to take effect at the time of application. Do not use if packaging is damaged or unintentionally opened prior to use.

STERILISATION**OptiLube** is supplied sterile.

Sterilised with gamma radiation after the packaging process.

Do not re-sterilise.

STORAGE AND HANDLING

The scale on the syringe is for the orientation of the user. It does not have a measurement function.

Store between 5-30°C (41-86°F) until expiry date.

Keep dry and out of direct sunlight.

Unopened packs have a 3-year shelf life. Do not use after expiry date.

DISPOSAL**OptiLube** must be disposed of according to local policies and procedures, including any remaining fluid in the syringe.

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