optilube.<sup>TM</sup>  
active CHG free

Instructions for use

## PRODUCT REF/CONTENTS

1167 **OptiLube Active CHG Free** 6ml Prefilled Syringe  
1168 **OptiLube Active CHG Free** 11ml Prefilled Syringe

This is an **OptiLube Active CHG Free** product information leaflet for the user. **OptiLube Active CHG Free** has been prescribed for you by a Healthcare Professional. Do not share it with anyone else. Read all of the product information leaflet prior to use and retain for your information. If you are unsure about how to use this product, or whether it is suitable for you, please seek medical advice.

## PRODUCT DESCRIPTION

**OptiLube Active CHG Free** (with **Lidocaine**) is a sterile, single-use, pre-filled syringe with water-soluble jelly available in 6ml and 11ml syringes.

## INGREDIENTS

100g of **OptiLube Active CHG Free** contains:  
Active Ingredient: • 2g **Lidocaine Hydrochloride**  
Other Ingredients: • Purified Water; • Propylene Glycol; • Hydroxyethyl cellulose.

## INTENDED PURPOSE

**OptiLube Active CHG Free** is a medical device that helps to prevent trauma being caused to the patient during catheterisation procedures or other urethral procedures by effective lubrication. Additionally, **OptiLube Active CHG Free** contains the following ancillary medicinal substance; local anaesthetic (**Lidocaine**) to help reduce pain for the patient.

## INDICATION FOR USE

**OptiLube Active CHG Free** is designed to lubricate the urethra prior to insertion of a urethral catheter and other urological medical devices including cystoscopes.

## CONTRAINDICATIONS

- OptiLube Active CHG Free** must not be used in patients:
- With known allergies or hypersensitivity to any of the ingredients.
  - Who have ever had a reaction to **Lidocaine**.
  - Who have damaged or bleeding mucous membranes because of the risk of systemic absorption of the **Lidocaine**.

Do not use in children below 2 years of age.  
Do not use for intravenous (IV) or intramuscular (IM) injections.  
Do not use orally. If swallowed, please seek medical advice.  
Do not use in eyes. If it comes into contact with the eyes, please seek medical advice.

## WARNINGS AND PRECAUTIONS

**OptiLube Active CHG Free** must only be used under the supervision of Healthcare Professionals in accordance with local guidelines, policies and procedures.

**OptiLube Active CHG Free** is not suitable if a patient:

- Has a damaged or bleeding urethra.
- Care should be taken when using the jelly with patients who:
- Have heart problems or taking medication for treating irregular heartbeat;
  - Have liver problems;
  - Are epileptic.

**OptiLube Active CHG Free** must not be used at the same time as other medical devices or medicines containing **Lidocaine** or other local anaesthetic.

Patients may experience a slight stinging on application of the jelly, but this usually stops once the anaesthetic starts to work. Encourage patients to report any reaction to the jelly.

Keep out of reach of children.

Ensure **OptiLube Active CHG Free** is suitable for the intended use and compatible with other medical devices to be used in conjunction with it.

This is a single-use device. Re-use of this device may result in patient infection/ cross-contamination.

Re/sterilisation, reprocessing, cleaning and disinfection may also compromise the product characteristics, resulting in trauma or infection to the patient.

**OptiLube Active CHG Free** does not contain natural rubber latex.

## PREGNANCY AND LACTATION

Only use during pregnancy and breastfeeding under the direction of a doctor. Always ask the patient if they are pregnant prior to using the jelly.

## EFFECTS ON THE ABILITY TO DRIVE AND USE MACHINERY

The ability to drive and operate machinery may be slightly impaired after the use of lubricant jelly with **Lidocaine**. If affected, patients should be advised not to drive or use machinery.

## INTERACTIONS OF MEDICINAL SUBSTANCES AND SUBSTANCES THAT ARE ABSORBED

Depending on the absorption of **Lidocaine**, these interactions can be seen when used with the following medications:

Propranolol: Reduction in plasma clearance of **Lidocaine**;

Cimetidine: Reduction in plasma clearance of **Lidocaine**;

Antiarrhythmic products: Increase in the toxicity of **Lidocaine**;

Phenytin or barbiturates: Reduction in plasma levels of **Lidocaine**.

Specified interactions can be seen in the long-term use and repeated high doses. When administered as recommended, there is no clinically significant interactions reported.

## UNDESIRABLE SIDE EFFECTS

Like any medication, **OptiLube Active CHG Free** may cause side effects in some people.

Side effects to this product must be documented in patient records.

In rare cases local hypersensitivity reactions may occur, such as redness, stinging, blistering or itching and/or systemic reactions to **Lidocaine**. There is also a risk of severe reactions including drop in blood pressure, dizziness, nausea, shortness of breath, bradycardia, convulsions and anaphylactic shock. If this occurs, stop using **OptiLube Active CHG Free** and consult your Healthcare Professional.

Any serious incident or malfunction which may affect the safety of **OptiLube Active CHG Free** should be reported to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

## DO dosage

Dosage recommendations:

FOR ADULTS: Max. 800 mg **Lidocaine** during a 24 hour period (x 3 syringe for 11 ml) (x 6 syringe for 6 ml),

FOR CHILDREN (between 2 and 15 years): Max. 0.3ml gel/kg B.W ( $\approx$ 6mg **Lidocaine**/kg) is recommended per procedure and no more than four doses should be administered within 24 hours.

In children below 2 years lubricant gel with **Lidocaine** must not be used.

## OVERDOSE

Excessive dosage or short intervals between doses can result in high plasma levels and serious adverse effects.

**OptiLube Active CHG Free** must not be used at the same time as any other medical device or medicines containing a local anaesthetic agent.

In the event of excessive absorption of **Lidocaine** into the bloodstream, symptoms may include central nervous system effects and cardiovascular reactions.

In the event of excessive use or use outside of these guidelines please seek medical advice.

## HOW TO USE

The decision as to which size (6ml or 11ml) of **OptiLube Active CHG Free** sterile lubricating jelly pre-filled syringe to be used is taken by the Healthcare Professional.

It is recommended that this product is used as part of an aseptic technique.

- Cleanse the opening to the urethra and surrounding area prior to use;
- Evenly peel back the paper backing and remove the sterile syringe;
- Remove the cap from the end of syringe;
- Apply a drop of jelly to urethral opening to make initial insertion easier;
- Insert the nozzle into the urethral opening and press the plunger of the syringe slowly to release the necessary amount of **OptiLube Active CHG Free**.

NOTE: **OptiLube Active CHG Free** should be applied to the urethra, not directly onto the device;

- The lubrication characteristics of the jelly start to take effect at the time of application. The onset of anaesthetic effect takes 3-5 minutes.

Do not use if packaging is damaged or unintentionally opened prior to use.

## STERILISATION

**OptiLube Active CHG Free** is supplied sterile.

Sterilised with gamma radiation after the packaging process.

Do not re-sterilise.

## STORAGE AND HANDLING

The scale on the syringe is for the orientation of the user. It does not have a measurement function.

Store between 5-30°C (41-86°F) until expiry date.

Keep dry and out of direct sunlight.

Unopened packs have a 3-year shelf life. Do not use after expiry date.

## DISPOSAL

**OptiLube Active CHG Free** must be disposed of according to local policies and procedures, including any remaining fluid in the syringe.

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1167 **OptiLube Active CHG Free** sans gluconate de chlorhexidine 6 ml Seringue préremplie  
1168 **OptiLube Active CHG Free** sans gluconate de chlorhexidine 11 ml Seringue préremplie

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Decomplex AG
Optimum Medical Solutions Ltd
Optimum Medical Solutions GmbH,
c/o Grant Thornton AG,