Black

BACK COVER

# optilube.<sup>TM</sup> active

REF 1160  
REF 1161

Rev. date: 11.12.2023

CE 2292

Does not contain natural rubber latex.  
No contiene látex de caucho natural.

STERILE R

MD

LATEX  
FREE

STERILIZE

PAPER ONLY

CH REP

MedEnvoy Switzerland  
Gotthardstrasse 28  
6302 Zug, Switzerland

Sponsor:

Optimum Medical Australasia Pty Ltd,  
Level 4, 470 Church Street, North Parramatta,  
New South Wales 2151, Australia

Istem Medikal Tıbbi Cihaz ve San. Tic. Ltd. Şti.,  
Anadolu OSB Mah. 29. Ekim Cad.,  
No: 41 Maltepe, Sincan, Ankara, Turkey

Made in Turkey

MedEnvoy Global B.V.  
Prinses Margrietplantsoen 33 - Suite 123  
2595 AM The Hague, The Netherlands

optimummedical®

UK Responsible Person and Distributor:  
Optimum Medical Solutions Ltd  
Tennant Hall, Blenheim Grove,  
Leeds, LS2 9ET, UK  
T: +44(0) 845 643 5479  
[enquiries@optimummedical.co.uk](mailto:enquiries@optimummedical.co.uk)  
[www.optimummedical.co.uk](http://www.optimummedical.co.uk)

1160-1161-A0140-Dec23-V3.0

Black

FRONT COVER

# optilube.<sup>TM</sup>

active

## Instructions for use

## GBR INSTRUCTIONS FOR USE

### PRODUCT REF/CONTENTS

1160 OptiLube Active 6ml Prefilled Syringe

1161 OptiLube Active 11ml Prefilled Syringe

This is an OptiLube Active product information leaflet for the user.

OptiLube Active has been prescribed for you by a Healthcare Professional. Do not share it with anyone else. Read all of the product information leaflet prior to use and retain for your information. If you are unsure about how to use this product, or whether it is suitable for you, please seek medical advice.

### PRODUCT DESCRIPTION

OptiLube Active (with Lidocaine and Chlorhexidine) is a sterile, single-use, pre-filled syringe with water-soluble jelly available in 6ml and 11ml syringes.

### INGREDIENTS

100g of OptiLube Active contains:

#### Active Ingredients

- 0.05g Chlorhexidine Gluconate; • 2g Lidocaine Hydrochloride.

#### Other Ingredients

- Purified Water; • Propylene Glycol; • Hydroxyethyl cellulose;
- Methyl Hydroxybenzoate; • Propyl Hydroxybenzoate.

### INTENDED PURPOSE

OptiLube Active is a medical device that helps to prevent trauma being caused to the patient during catheterisation procedures or other urethral procedures by effective lubrication. Additionally, OptiLube Active contains the following ancillary medicinal substances; antiseptic (Chlorhexidine) and local anaesthetic (Lidocaine) to help reduce the risk of infection and pain for the patient.

### INDICATION FOR USE

OptiLube Active is designed to lubricate the urethra prior to insertion of a urethral catheter and other urological medical devices including cystoscopes.

### CONTRAINDICATIONS

OptiLube Active must not be used in patients:

- With known allergies or hypersensitivity to any of the ingredients.
- Who have ever had a reaction to Lidocaine or Chlorhexidine.
- Who have damaged or bleeding mucous membranes because of the risk of systemic absorption of the Lidocaine and Chlorhexidine.

Do not use in children below 2 years of age.

Do not use for intravenous (IV) or intramuscular (IM) injections.

Do not use orally. If swallowed, please seek medical advice.

Do not use in eyes. If it comes into contact with the eyes, please seek medical advice.

### WARNINGS AND PRECAUTIONS

OptiLube Active must only be used under the supervision of Healthcare Professionals in accordance with local guidelines, policies and procedures.

OptiLube Active is not suitable if a patient:

- Has a damaged or bleeding urethra.

Care should be taken when using the jelly with patients who:

- Have heart problems or taking medication for treating irregular heartbeat;
- Have liver problems;
- Are epileptic.

**OptiLube Active** must not be used at the same time as other medical devices or medicines containing **Lidocaine** or other local anaesthetic.

Patients may experience a slight stinging on application of the jelly, but this usually stops once the anaesthetic starts to work. Encourage patients to report any reaction to the jelly.

Keep out of reach of children.

Ensure **OptiLube Active** is suitable for the intended use and compatible with other medical devices to be used in conjunction with it.

This is a single-use device. Re-use of this device may result in patient infection/cross-contamination.

Re/sterilisation, reprocessing, cleaning and disinfection may also compromise the product characteristics, resulting in trauma or infection to the patient.

**OptiLube Active** does not contain natural rubber latex.

#### **PREGNANCY AND LACTATION**

Only use during pregnancy and breastfeeding under the direction of a doctor.

Always ask the patient if they are pregnant prior to using the jelly.

#### **EFFECTS ON THE ABILITY TO DRIVE AND USE MACHINERY**

The ability to drive and operate machinery may be slightly impaired after the use of lubricant jelly with **Lidocaine**. If affected, patients should be advised not to drive or use machinery.

#### **INTERACTIONS OF MEDICINAL SUBSTANCES AND SUBSTANCES THAT ARE ABSORBED**

Depending on the absorption of **Lidocaine**, these interactions can be seen when used with the following medications:

Propranolol: Reduction in plasma clearance of **Lidocaine**;

Cimetidine: Reduction in plasma clearance of **Lidocaine**;

Antiarrhythmic products: Increase in the toxicity of **Lidocaine**;

Phenytoin or barbiturates: Reduction in plasma levels of **Lidocaine**.

Specified interactions can be seen in the long-term use and repeated high doses.

When administered as recommended, there is no clinically significant interactions reported.

#### **UNDESIRABLE SIDE EFFECTS**

Like any medication, **OptiLube Active** may cause side effects in some people. Side effects to this product must be documented in patient records.

In rare cases local hypersensitivity reactions may occur, such as redness, stinging, blistering or itching and/or systemic reactions to **Lidocaine** and/or **Chlorhexidine**.

There is also a risk of severe reactions including drop in blood pressure, dizziness, nausea, shortness of breath, bradycardia, convulsions and anaphylactic shock. If this occurs, stop using **OptiLube Active** and consult your Healthcare Professional.

Any serious incident or malfunction which may affect the safety of **OptiLube Active** should be reported to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

## **DOSAGE**

Dosage recommendations:

**FOR ADULTS:** Max. 800 mg **Lidocaine** during a 24 hour period (x 3 syringe for 11 ml) (x 6 syringe for 6 ml),

**FOR CHILDREN** (between 2 and 15 years): max. 0.3ml gel/kg B.W ( $\cong$ 6mg **Lidocaine**/kg) is recommended per procedure and no more than four doses should be administered within 24 hours.

In children below 2 years lubricant gel with **Lidocaine** must not be used.

## **OVERDOSE**

Excessive dosage or short intervals between doses can result in high plasma levels and serious adverse effects.

**OptiLube Active** must not be used at the same time as any other medical device or medicines containing a local anaesthetic agent.

In the event of excessive absorption of **Lidocaine** into the bloodstream, symptoms may include central nervous system effects and cardiovascular reactions.

In the event of excessive use or use outside of these guidelines please seek medical advice.

## **HOW TO USE**

The decision as to which size (6ml or 11ml) of **OptiLube Active** sterile lubricating jelly pre-filled syringe to be used is taken by the Healthcare Professional. It is recommended that this product is used as part of an aseptic technique.

- Cleanse the opening to the urethra and surrounding area prior to use;
- Evenly peel back the paper backing and remove the sterile syringe;
- Remove the cap from the end of syringe;
- Apply a drop of jelly to urethral opening to make initial insertion easier;
- Insert the nozzle into the urethral opening and press the plunger of the syringe slowly to release the necessary amount of **OptiLube Active**.

**NOTE:** **OptiLube Active** should be applied to the urethra, not directly onto the device;  
• The lubrication characteristics of the jelly start to take effect at the time of application.

The onset of anaesthetic effect takes 3-5 minutes.

Do not use if packaging is damaged or unintentionally opened prior to use.

## **STERILISATION**

**OptiLube Active** is supplied sterile.

Sterilised with gamma radiation after the packaging process.

Do not resterilise.

## **STORAGE AND HANDLING**

The scale on the syringe is for the orientation of the user.

It does not have a measurement function.

Store between 5-30°C (41-86°F) until expiry date.

Keep dry and out of direct sunlight.

Unopened packs have a 3-year shelf life. Do not use after expiry date.

## DISPOSAL

OptiLube Active must be disposed of according to local policies and procedures, including any remaining fluid in the syringe.

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