



constantly creating better

THE UgoBOX™

a fully equipped catheterisation kit for the insertion and removal of urethral catheters



We're specialists in urinary drainage and fixation products.

Find out more about our Ugo Urology range by visiting www.optimummedical.co.uk.



Optimum Medical Solutions Ltd
Tennant Hall, Blenheim Grove,
Leeds, LS2 9ET, UK
T: +44 (0) 845 643 5479
enquiries@optimummedical.co.uk
www.optimummedical.co.uk

EC REP

MT Promed Consulting GmbH
Altenhofstrasse 80
66386 St. Ingbert Germany



OMS-UB-3000-May20



THE **ugobox**.™



The **Ugo Box** is a fully equipped catheterisation kit for the insertion and removal of urethral catheters. Each box comes equipped with a choice of **Ugo Foley Catheter**, **Ugo Leg Bag** and **OptiLube** catheterisation jelly syringe.

The **Ugo Box** has been designed to support best practice aseptic catheterisation, promote efficiency, and provide a cost-effective solution in one handy pack.

The contents of the **Ugo Box** must be used in accordance with local policies and procedures and best practice guidance by suitably qualified and competent healthcare practitioners.

THE UGO BOX CONTAINS:

- 1 1X **Ugo Pro** (catheterisation procedure pack with 20ml saline and 10ml empty syringe)
- 2 1X **Ugo Foley Catheter** (all silicone)
- 3 1X **OptiPure** (10ml sterile water syringe)
- 4 1X 11 ml catheterisation jelly syringe (**OptiLube** or **OptiLube Active CHG Free**)
- 5 1X **Ugo Leg Bag** (1C – short tube, 500ml, lever tap or 2C – long tube, 500ml, lever tap)
- 6 1X **Ugo Fix Gentle** (catheter clip)
- 7 1X pair **Ugo Fix Leg Bag Straps**



ORDERING INFORMATION

We understand the importance of having access to a choice. This is why we provide the **Ugo Box** in 12 variations.

PR Code	PIP Code	Ugo Foley Catheter FR size	Ugo Leg Bag	11ml catheterisation jelly syringe
3450	415-0660	12FR	Ugo 1C	OLA CHG Free*
3451	415-0678	14FR	Ugo 1C	OLA CHG Free
3452	415-0686	16FR	Ugo 1C	OLA CHG Free
3453	415-0694	12FR	Ugo 1C	OptiLube
3454	415-0702	14FR	Ugo 1C	OptiLube
3455	415-0710	16FR	Ugo 1C	OptiLube
3456	415-0728	12FR	Ugo 2C	OLA CHG Free
3457	415-0736	14FR	Ugo 2C	OLA CHG Free
3458	415-0744	16FR	Ugo 2C	OLA CHG Free
3459	415-0769	12FR	Ugo 2C	OptiLube
3460	415-0751	14FR	Ugo 2C	OptiLube
3461	415-0777	16FR	Ugo 2C	OptiLube

Please note: The **Ugo Foley Catheter** balloon size is 10ml.

* **OptiLube Active CHG Free**

GIVING CLINICIANS A CHOICE OF CATHETERISATION GEL

OptiLube

OptiLube is a sterile, water-soluble lubricating jelly without active ingredients. The lack of active ingredients makes it ideal for catheterising patients with allergies or contraindications.

OptiLube Active CHG Free

OptiLube Active CHG Free is a sterile, water-soluble lubricating jelly with **Lidocaine 2%**. The addition of **Lidocaine 2%** helps to reduce the risk of infection.

KEY FEATURES

- Cost effective.
- Recyclable packaging.
- Reduces waste - by providing all the items required in one box, this reduces stockpiling, over-prescribing and excessive packaging.
- Promotes efficiency - reduces time spent setting up the procedure, resulting in more time to spend with the patient.
- The fully comprehensive IFU supports both HCP and patient education of the procedure and the components included.
- Compliant with National Guidelines and evidence-based practices.
- Adheres to all local policies and procedures or techniques.
- 100% latex-free.

UGO LEG BAG 1C OR 2C?

The **Ugo Box** has 2 variations of **Ugo Leg Bags** to choose from:

Ugo 1C (500ml, 5cm short tube, lever tap) or **Ugo 2C** (500ml, 25cm long tube, lever tap).

Both **Ugo 1C** and **Ugo 2C** can remain attached to a patient's catheter for up to 7 days. Each **Ugo Leg Bag** is sterile pouch packaged and manufactured from a durable yet soft medical-grade plastic material for strength and reduced noise when moving.