



optiPure™

PRE-FILLED STERILE WATER SYRINGES


Pre-filled syringes containing sterile water are typically used to inflate the catheter balloon of a Foley indwelling catheter.

OptiPure is a class I sterile medical device. Available in two formats for added patient comfort; sterile water pre-filled syringe and sterile water pre-filled syringe with 10% glycerine, designed to maintain a higher level of catheter balloon inflation.

OptiPure features and benefits:

- Use of syringe format allows for a constant flow control during application
- Measuring facility on pre-filled syringe ensures an accurate application of fluid
- Sterilised inside and out, with easy peel packaging; enabling safe use in sterile areas
- Available in two formats for added patient comfort.

Code	Description	Inner Pack	Outer Pack
1128	10ml Sterile Water Pre-filled Syringe	Box of 20	Case of 200
1131	10ml Sterile Water Pre-filled Syringe with 10% Glycerine	Box of 20	Case of 200



LATEX FREE

STERILE R CE 0473



optimummedical®

Unit 5, First Floor, Armley Court,
Leeds, LS12 2LB, UK
Tel.: +44(0)845 643 5479,
Fax: +44(0)114 238 3826
enquiries@optimummedical.co.uk
optimummedical.co.uk



Company Registration No: 6779883