


constantly creating better

Ugo™
Foley catheter



We're specialists in
healthcare lubrication and
urinary drainage products.

Find out more about our Ugo Urology and
OptiLube catheterisation jelly products by visiting
www.optimummedical.co.uk.



Optimum Medical Solutions Ltd
Tennant Hall, Blenheim Grove,
Leeds, LS2 9ET, United Kingdom
T: +44(0) 845 643 5479
enquiries@optimummedical.co.uk
www.optimummedical.co.uk

EC REP

MT Promed Consulting GmbH
Altenhofstrasse 80
66386 St. Ingbert Germany



OMS-UFC-200-Jan20



The **Ugo Foley Catheter** is designed for urethral catheterisation of the urinary bladder. When fitted by a qualified healthcare professional following a clinical assessment of a patient, it may be used for continuous bladder drainage, instillation of appropriate therapeutic agents into the bladder, or bladder irrigation.

Patients requiring long-term catheters should be continuously assessed and have them removed when not required.

KEY FEATURES

- Made from 100% silicone.
- Latex-free.
- Integrated balloon.
- Standard length (420mm).
- Up to 12 weeks use.
- X-ray detectable line.
- Includes 2 catheter product record stickers to support best practice.



IMPORTANT:

- Catheters should always be inserted using an aseptic technique.
- Catheters should be secured using a fixation device, such as the **Ugo Fix Gentle** (catheter clip) or **Ugo Fix Catheter Strap**.

Find out more at www.optimummedical.co.uk.

ORDERING INFORMATION

Each **Ugo Foley Catheter** includes:

- 1x **OptiPure** 10ml sterile water syringe
- 1x 10ml **Ugo** empty syringe, and full instructions for use.

PR Code	PIP Codes	Product Name	Size
UFC12	414-3558	Ugo Foley Catheter	12FR
UFC14	414-3541	Ugo Foley Catheter	14FR
UFC16	414-3533	Ugo Foley Catheter	16FR
UFC18	414-3525	Ugo Foley Catheter	18FR

Male and female catheterisation animations and a supporting booklet are available, request a free DVD and booklet by emailing marketing@optimummedical.co.uk.

CATHETERISATION JELLY

To give clinicians more choice, we offer a range of lubricating gels for catheterisation procedures, all provided in a familiar syringe format.



optilube™ active with Lidocaine 2% and Chlorhexidine 0.05%.	PR Code	PIP Codes	Description	Inner Pack	Outer Pack
	1160	366-2061	OptiLube Active 6ml syringe	Box of 10	Case of 200
	1161	366-2079	OptiLube Active 11ml syringe	Box of 10	Case of 200
optilube™ active CHG free with Lidocaine 2%.	PR Code	PIP codes	Description	Inner Pack	Outer Pack
	1167	404-0663	OptiLube Active CHG Free 6ml syringe	Box of 10	Case of 200
	1168	404-0671	OptiLube Active CHG Free 11ml syringe	Box of 10	Case of 200
optilube™ without active ingredients.	PR Code	PIP codes	Description	Inner Pack	Outer Pack
	1125	373-4357	OptiLube 6ml syringe	Box of 25	Case of 250
	1126	373-4365	OptiLube 11ml syringe	Box of 25	Case of 250