

FFP3

NEW!

Vertical fold-flat disposable face mask with exhalation valve



- ✓ Compliant with guidelines EN 149:2001 & A1:2009.
- ✓ Exhalation valve reduces moisture and heat build-up.
- ✓ Cushion-fit lining offers increased comfort and security when worn.
- ✓ Four-point headband attachment for extra security.

PR CODE	PRODUCT NAME	INNER PACK	OUTER PACK
HY8232	FFP3 mask	Box of 20	Case of 20

HY® | Manufactured by Handanhy for Optimum Medical

To find out more information, email: enquiries@optimummedical.co.uk

Or give us a call directly on: **+44 (0) 113 263 3849**



optimummedical® | constantly creating better

Tennant Hall, Blenheim Grove, Leeds, LS2 9ET, UK