



constantly creating better

Ugo™
leg bag



We're specialists in
urology products.

Find out more about our Ugo Urology and
OptiLube catheterisation jelly products by visiting
www.optimummedical.co.uk.



Optimum Medical Solutions Ltd
Tennant Hall, Blenheim Grove,
Leeds, LS2 9ET, United Kingdom
T: +44 (0) 845 643 5479
enquiries@optimummedical.co.uk
www.optimummedical.co.uk

EC REP

MT Promed Consulting GmbH
Altenhofstrasse 80
66386 St. Ingbert Germany



OMS-UDB-105-Aug20





- All **Ugo Leg Bags** are sterile pouch packaged which is printed with key user information in a clear and simple format.
- All **Ugo Leg Bags** can remain attached to an indwelling catheter for up to 7 days.

At **Optimum Medical** we believe we have a responsibility to provide enough information that those using our products have a safe, effective and positive experience*.

Ugo Leg Bags are designed with a range of features to meet individual patient and clinical requirements**, and each is provided with a comprehensive patient and carer user guide to make the experience using our products the best it can be.

* Darzi's Definition of High Quality Care
 ** In line with NICE (2012) and epic3 best practice guidelines

KEY STRENGTHS

- Transparent universal stepped connector for secure fitting into all indwelling catheters.
- Needle free sprung sample port which locks syringe for safe sample collection.
- Non return valve to prevent urine returning up the drainage tube.
- 4 strengthened eyelets to ensure level fixation using enclosed **Ugo Fix Leg Bag Straps**.
- Discreet bag printing with unique naming system for quick and easy re-ordering.
- Purpose measured soft silicone outlet tubing for perfect kink-free link drainage connection.
- Natural leg shape design
- Choice of fabric or non fabric backing depending on clinical environment and patient choice
- Durable but soft medical grade plastic material for strength and reduced noise when moving



AVAILABLE IN:

Lever tap

or

T tap

Our **Ugo Pro Leg Bag** range are available with a pair of disposable gloves and apron in each pouch to reduce infection risk to both patients and their carers. (Please see ordering information on the next page)

FIXATION DEVICES

To secure and support catheter and drainage systems



Find out more:
www.optimummedical.co.uk

The range includes:

- Ugo Fix Gentle** (catheter clip)
- Ugo Fix Catheter Strap** (catheter fixation strap)
- Ugo Fix Sleeve** (leg bag holder)
- Ugo Fix Bag Belt** (leg bag holder)
- MultiFix Gentle** (universal securement device for tubing and urinary catheters)

ORDERING INFORMATION

Ugo Leg Bags (Community) - Fabric backing – 1 pair of **Ugo Fix Leg Bag Straps** in every box

PR CODE	PR NAME	PIP CODE	BAG VOL	TUBE LENGTH	TAP DESIGN	INNER PACK	OUTER CASE
1C/500/ST/L	Ugo 1C	400-8546	500ml	Short Tube	Lever Tap	Box of 10	8 boxes of 10
2C/500/LT/L	Ugo 2C	400-8603	500ml	Long Tube	Lever Tap	Box of 10	8 boxes of 10
3C/500/ST/T	Ugo 3C	400-8611	500ml	Short Tube	T tap	Box of 10	8 boxes of 10
4C/500/LT/T	Ugo 4C	400-8629	500ml	Long Tube	T tap	Box of 10	8 boxes of 10
5C/350/DI/L	Ugo 5C	400-8637	350ml	Direct Inlet	Lever Tap	Box of 10	8 boxes of 10
6C/350/ST/L	Ugo 6C	400-8645	350ml	Short Tube	Lever Tap	Box of 10	8 boxes of 10
7C/750/ST/L	Ugo 7C	400-8652	750ml	Short Tube	Lever Tap	Box of 10	8 boxes of 10
8C/750/LT/L	Ugo 8C	400-8660	750ml	Long Tube	Lever Tap	Box of 10	8 boxes of 10

Ugo Leg Bags (Hospital) - No fabric backing – 1 pair of **Ugo Fix Leg Bag Straps** in every pouch

PR CODE	PR NAME	BAG VOL	TUBE LENGTH	TAP DESIGN	INNER PACK	OUTER CASE
1H/500/ST/L	Ugo 1H	500ml	Short Tube	Lever Tap	Box of 10	8 boxes of 10
2H/500/LT/L	Ugo 2H	500ml	Long Tube	Lever Tap	Box of 10	8 boxes of 10
3H/500/ST/T	Ugo 3H	500ml	Short Tube	T tap	Box of 10	8 boxes of 10
4H/500/LT/T	Ugo 4H	500ml	Long Tube	T tap	Box of 10	8 boxes of 10

Ugo Fix Leg Bag Straps

PR CODE	PIP CODE	DESCRIPTION	INNER PACK	OUTER CASE
3011	400-8538	Ugo Fix Leg Bag Straps	10 pairs *	Case of 200

* Each pair contains one 75cm and one 45cm leg bag strap

Ugo Pro Leg Bags –Fabric backing –includes disposable gloves and apron, and 1 pair of **Ugo Fix Leg Bag Straps** in every box

PR CODE	PR NAME	PIP CODE	BAG VOL	TUBE LENGTH	TAP DESIGN	INNER PACK	OUTER CASE
3410	Ugo Pro 1C	415-9554	500ml	Short Tube	Lever Tap	Box of 10	8 boxes of 10
3411	Ugo Pro 2C	415-9547	500ml	Long Tube	Lever Tap	Box of 10	8 boxes of 10
3412	Ugo Pro 3C	415-9570	500ml	Short Tube	T tap	Box of 10	8 boxes of 10
3413	Ugo Pro 4C	415-9562	500ml	Long Tube	T tap	Box of 10	8 boxes of 10



References:

- Quality governance in the NHS - a guide for provider boards (2011) https://www.gov.uk/government/uploads/system/uploads/attachment_data/file/216321/dh_125239.pdf
- National Institute for Health and Care Excellence (2012) Guidance 139. Infection: Prevention and Control of Healthcare-associated Infections in Primary and Community Care. London: National Institute for Health and Care Excellence; Available at HYPERLINK "http://www.nice.org.uk" www.nice.org.uk.
- Loveday HP, Wilson JA, Pratt RJ et al. (2013) Epic3: National Evidence-Based guidelines for Preventing Healthcare-Associated Infections in NHS Hospitals in England. Journal of Hospital Infection; 86S1:51-570.